



# pK<sub>a</sub> measurement using 50 - 100µL aliquots from DMSO stock, no chromophore required

Karl Box<sup>1</sup>, John Comer<sup>1</sup>, Tom Gravestock<sup>1</sup>, Jon Mole<sup>2</sup>

<sup>1</sup>Sirius Analytical, Forest Row, East Sussex, RH18 5DW, UK.  
<sup>2</sup>Sirius Analytical Inc. Suite A, 1825 Swarthmore Avenue, Lakewood, NJ 08701

## Abstract

**Purpose:** pH-metric pK<sub>a</sub> measurement requires a weighed sample of pure drug. For titrations done in 1mL of solution, the amount required is so small that it is difficult to weigh accurately. A method was therefore sought for taking samples as aliquots of solution in DMSO stock.


**Methods:** Potentiometric titration is a powerful method for pK<sub>a</sub> measurement because it is based only on pH measurement (hence the name pH-metric). It can therefore measure pK<sub>a</sub>(s) of all drugs, including those such as aliphatic amines or carboxylic acids that show no significant UV absorbance associated with the ionizable group. Though powerful, the method has been difficult to apply for drug samples from discovery because the weight required for titrating in 10mL volume (2 – 5 mg) may be too high. A new instrument has been developed for titrating in 1 mL volume of solution, making it possible to do pH-metric titrations using 0.2mg of sample. Because such low amounts are difficult to weigh, samples were prepared using 50-100µL aliquots of 10mM sample stock solution in DMSO.

**Results:** The goal of this study was to investigate the effect of volumes of DMSO introduced on the aqueous pK<sub>a</sub> value. 100µL of DMSO represents 10% of the volume of a titration in 1mL, and it is known that the presence of DMSO changes the apparent pK<sub>a</sub> (see pramoxine example in Figure 1a). This work shows that the average shift in pK<sub>a</sub> caused by 10% DMSO varies between 0.03 and 0.26, and that the shift is similar when pK<sub>a</sub>(s) are measured in aqueous solution (containing 10% DMSO), or when they are determined by extrapolation from successive measurements made in different ratios of water-methanol (also containing 10% DMSO). These shifts vary little from sample to sample, and could either be ignored, or applied to the measured result as a correction factor.

**Conclusion:** pK<sub>a</sub>(s) can be measured by the pH-metric method using 50-100µL aliquots of 10mM DMSO sample stocks in 1mL volumes. The shifts in pK<sub>a</sub> caused by the DMSO have been quantified.

## Experimental

Six drugs were selected for study. The ionizable groups of four of these (imipramine, captopril, gabapentin and pramoxine) were remote from UV-absorbing chromophores, typical of drugs whose pK<sub>a</sub> values cannot be measured by pH-UV methods. Samples were measured by the pH-metric method, under a variety of conditions shown in Table 1. A pK<sub>a</sub> result is obtained from titration data collected in a water/solvent mixture. All experiments were done at 25±1°C, in solutions with ionic strength adjusted to 0.15M using KCl. As well as pK<sub>a</sub> measured in water/methanol, pK<sub>a</sub> was measured in 0.15M aqueous KCl (with and without 10% DMSO) for the water-soluble compounds (pramoxine, captopril, gabapentin and famotidine).



Sample	pK <sub>a</sub> (1)	pK <sub>a</sub> (2)	Mean methanol Wt % per pair of points
Pramoxine: base with one pK <sub>a</sub>			
Extrapolated pK <sub>a</sub> , no DMSO	7.323		0 12 21 30 40 50
Extrapolated pK <sub>a</sub> , 10% DMSO	6.906		
ΔpK <sub>a</sub>	-0.217		
Imipramine: base with one pK <sub>a</sub>			11 21 30 40 50 60
Extrapolated pK <sub>a</sub> , no DMSO	9.467		
Extrapolated pK <sub>a</sub> , 10% DMSO	9.331		
ΔpK <sub>a</sub>	-0.334		
Flumequine: acid with one pK <sub>a</sub>			12 21 30 40 50 60
Extrapolated pK <sub>a</sub> , no DMSO	6.307		
Extrapolated pK <sub>a</sub> , 10% DMSO	6.375		
ΔpK <sub>a</sub>	+0.068		
Captopril: acid with two pK <sub>a</sub> s			0 12 21 31 40 50 60
Extrapolated pK <sub>a</sub> , no DMSO	9.807	3.485	
Extrapolated pK <sub>a</sub> , 10% DMSO	9.773	3.615	
ΔpK <sub>a</sub>	-0.034	+0.13	
Gabapentin: zwitterion with two pK <sub>a</sub> s			0 12 21 30 40 49 60
Extrapolated pK <sub>a</sub> , no DMSO	10.647	9.728	
Extrapolated pK <sub>a</sub> , 10% DMSO	10.84	3.82	
ΔpK <sub>a</sub>	-0.207	+0.092	
Famotidine: ampholyte with two pK <sub>a</sub> s			0 12 21 30 38 48 60
Extrapolated pK <sub>a</sub> , no DMSO	11.279	6.777	
Extrapolated pK <sub>a</sub> , 10% DMSO	11.162	6.512	
ΔpK <sub>a</sub>	-0.117	-0.265	

Table 1 Compounds studied and extrapolated pK<sub>a</sub> results.

The mean of all extrapolated acidic ΔpK<sub>a</sub> values is +0.03 ± 0.10

The mean of all extrapolated basic ΔpK<sub>a</sub> values is -0.26 ± 0.06

**Results** The aqueous pK<sub>a</sub> values of poorly water-soluble drugs are determined by Yasuda-Shedlovsky (Y-S) extrapolation from pK<sub>a</sub> values measured in solutions with different ratios of water to solvent (e.g. water/methanol). In these Y-S extrapolations, pK<sub>a</sub> + log[H<sub>2</sub>O] is plotted against 1/ε, where ε is the dielectric constant of the water/solvent mixture. The aqueous pK<sub>a</sub> is obtained by subtracting log[H<sub>2</sub>O] of water (1.755) from the y-axis intercept at 1/ε for water, as illustrated in figures 1a and 1b.

Aqueous pK<sub>a</sub> values determined by Yasuda-Shedlovsky extrapolation with and without DMSO present are shown in Table 1.

For each pair of points in Figure 1, the pK<sub>a</sub> (or pK<sub>a</sub>) measured without DMSO was subtracted from the equivalent value measured in the presence of 10% DMSO, and recorded as ΔpK<sub>a</sub>, as plotted vs. mean methanol content in figure 2. Note that ΔpK<sub>a</sub> values are negative for all basic groups, and smaller and more variable for acidic groups.




Figure 1 Six graphs showing Yasuda-Shedlovsky plots of data collected with and without 10% DMSO. pK<sub>a</sub> values for acidic groups are denoted in red, pK<sub>a</sub> values for basic groups are denoted in blue. Methanol content increases from left to right on the x-axis by approximately 10% between each pair of points.




Figure 2 ΔpK<sub>a</sub> per pair of measured pK<sub>a</sub> points.

ΔpK<sub>a</sub> values for acidic groups are denoted in red. ΔpK<sub>a</sub> values for basic groups are denoted in blue.

The mean of all acidic ΔpK<sub>a</sub> values is +0.15 ± 0.15

The mean of all basic ΔpK<sub>a</sub> values is -0.23 ± 0.08

**Conclusions** The goal of this study was to quantify the differences between pK<sub>a</sub> and pK<sub>a</sub> values measured with and without 10% DMSO present. From the extrapolated pK<sub>a</sub> values reported in Table 1 and the ΔpK<sub>a</sub> values in Figure 2, we draw the following conclusions:

Average ΔpK <sub>a</sub> for all acidic pK <sub>a</sub> values	0.15 ± 0.15	Average ΔpK <sub>a</sub> for all basic pK <sub>a</sub> values	-0.23 ± 0.08
Average ΔpK <sub>a</sub> for extrapolated acidic pK <sub>a</sub> values	0.03 ± 0.10	Average ΔpK <sub>a</sub> for extrapolated basic pK <sub>a</sub> values	-0.26 ± 0.06

These differences caused by 10% DMSO are not very large, and may be within acceptable limits for many applications. Alternatively, a correction factor could be applied as above. The practical implications are that pH-metric pK<sub>a</sub> measurements can be run in 1mL volumes using 100µL aliquots of 10mM sample in DMSO with the Sirius T3 instrument.

Sirius - Instruments and PhysChem Analytical Services

www.sirius-analytical.com

## Instrumentation

pK<sub>a</sub>, logP & logD,  
Solubility & dissolution  
Surface Activity Profiling  
Phospholipidosis

## Analytical Services

PhysChem properties – pK<sub>a</sub>, logP, logD, solubility, dissolution  
Solid state assays – XRPD, DSC, TGA, Raman  
Surface Tension – CMC, TSA, K<sub>AW</sub>

For more information please contact  
sirius@sirius-analytical.com

Sirius Analytical Ltd.  
Riverside, Forest Row Business Park,  
Forest Row, East Sussex, RH18 5DW, UK  
Phone: +44 1342 820720 Fax: +44 1342 820725  
Web: [www.sirius-analytical.com](http://www.sirius-analytical.com)